

## ReVel<sup>®</sup> Job Aid: Breath Mode and Type

**Table 1: Matrix of Active Controls in All Breath Mode and Breath Type Settings**

<i>Breath Type &amp; Mode</i>	<i>A/C + Volume</i>	<i>A/C + Pressure</i>	<i>A/C + PRVC</i>	<i>SIMV + Volume</i>	<i>SIMV + Pressure</i>	<i>SIMV + PRVC</i>	<i>CPAP + PS + Volume</i>	<i>CPAP + PS + Pressure</i>	<i>CPAP + PS + PRVC</i>	<i>A/C + Pressure</i>	<i>CPAP + PS Pressure</i>
<b>CONTROLS</b>											
<b>Bias Flow<sup>1</sup></b>	*	*	*	*	*	*	*	*	*	*	*
<b>Breath Rate</b>	*	*	*	*	*	*				*	
<b>Flow Term<sup>1</sup></b>		* <sup>2</sup>	* <sup>2</sup>	*	*	*	*	*	*	* <sup>2</sup>	*
<b>Insp. Time</b>	*	*	*	*	*	*			*	*	
<b>Leak Comp<sup>1</sup></b>	*	*	*	*	*	*	*	*	*	*	*

<sup>1</sup> ReVel<sup>®</sup>, CareFusion, and the CareFusion logo are trademarks or registered trademarks of CareFusion or one of its subsidiaries. All other trademarks are the property of their respective owners.

<i>Breath Type &amp; Mode</i>	<i>A/C + Volume</i>	<i>A/C + Pressure</i>	<i>A/C + PRVC</i>	<i>SIMV + Volume</i>	<i>SIMV + Pressure</i>	<i>SIMV + PRVC</i>	<i>CPAP + PS + Volume</i>	<i>CPAP + PS + Pressure</i>	<i>CPAP + PS + PRVC</i>	<i>A/C + Pressure</i>	<i>CPAP + PS Pressure</i>
<b>CONTROLS</b>											
<b>O<sub>2</sub></b>	*	*	*	*	*	*	*	*	*	*	*
<b>PC Flow Term<sup>1</sup></b>		*	*		*	*				*	
<b>PEEP</b>	*	*	*	*	*	*	*	*	*	*	*
<b>Pres Trigger<sup>1</sup></b>	*	*	*	*	*	*	*	*	*	*	*
<b>Pres. Control</b>		*			*					*	
<b>Pres. Support</b>				*	*	*	*	*			*
<b>Rise Time<sup>1</sup></b>		*	*	*	*	*	*	*	*	*	*

2 ReVel<sup>®</sup>, CareFusion, and the CareFusion logo are trademarks or registered trademarks of CareFusion or one of its subsidiaries. All other trademarks are the property of their respective owners.

<i>Breath Type &amp; Mode</i>	<i>A/C + Volume</i>	<i>A/C + Pressure</i>	<i>A/C + PRVC</i>	<i>SIMV + Volume</i>	<i>SIMV + Pressure</i>	<i>SIMV + PRVC</i>	<i>CPAP + PS + Volume</i>	<i>CPAP + PS + Pressure</i>	<i>CPAP + PS + PRVC</i>	<i>A/C + Pressure</i>	<i>CPAP + PS Pressure</i>
<b>CONTROLS</b>											
<b>Sensitivity</b>	*	*	*	*	*	*	*	*	*	*	*
<b>Tidal Volume</b>	*		*	*		*			*		
<b>Time Term<sup>1</sup></b>				*	*	*	*	*			*
	<b>Intubated Patient Interface</b>									<b>Non-Invasive Patient Interface (NPPV)</b>	

<sup>1</sup> Control accessed while in normal ventilation modes through the Extended Features menus.

<sup>2</sup> Available when Pressure Control Flow Termination (PC Flow Term) is on.

**Caution:** Federal Law restricts this device to sale, distribution and use by or on the order of a physician.